


# SDMA

Research Diagnostics

**An introduction**





# Table of contents

What is SDMA?	3
Features & Benefits	4
Testing Methodologies	5
Kidneys in Action	6
Creatinine Connection	7
LCMS — The Gold Standard	8
Human Clinical Studies; partnering with Yale University	9
Resources	10

# What is SDMA?

## Serum symmetric dimethylarginine (SDMA)

SDMA stands for symmetric dimethylarginine. SDMA is an amino acid produced in the body when protein is broken down and then excreted through the kidneys. It is measured to detect kidney disease.

- SDMA is a methylated form of arginine in intracellular proteins of all nucleated cells
- SDMA is released into circulation when intracellular proteins are processed
- SDMA is excreted by the kidneys
- Stable production of SDMA is part of daily cell activity
- It is a renal biomarker shown to correlate with glomerular filtration rate (GFR)

## An earlier indicator of renal failure

Historically, estimation of renal function has depended on blood urea nitrogen (BUN) and serum creatinine. Creatinine is a crude estimate of renal function as it doesn't elevate until significant (~75%) loss of function occurs. **SDMA identifies at only 40% loss of kidney function and as early as 25% loss.\*** Compare this to creatinine testing which detects once 75% of kidney deterioration has occurred and is an earlier indicator of renal failure than creatinine/BUN.

# Features and benefits

- SDMA is a marker for kidney function
- SDMA is earlier, more sensitive, and more reliable than creatinine
- SDMA correlates with GFR
- SDMA is a sensitive indicator of kidney function that **detects as little as 25% loss of function.**
- SDMA is more reliable than creatinine as an indicator of kidney function because it is not influenced by common confounding conditions.
- Compared to BUN, SDMA is a more reliable and sensitive indicator of kidney function in animals. BUN can be influenced by decreased production of urea in liver disease or anorexia and by increased production with high-protein meals or gastrointestinal bleeding. This contrasts with SDMA, which changes only with changes in GFR.
- SDMA is an earlier indicator of progressive kidney function loss, often increasing before other parameters.
- An increased SDMA may also serve as an indicator of concurrent diseases that may have a secondary impact on kidney function.

Benefit	Mouse	Rat	Dog	NHP
Biomarker for kidney function	•	•	•	
Correlates with GFR	•	•	•	
Earlier indicator than Creatinine			•	
More sensitive than Creatinine			•	
Not subject to extrarenal factors	•	•	•	



# Testing methodologies

## Validation

- **SDMA has been validated in rats, and used in mice (basic science)** — LCMS
- **SDMA is validated in dogs, cats, and horses (clinical studies)**— LCMS, Beckman (Immunoassay), In house diagnostics (Immunoassay - Catalyst ®)
- **SDMA is highly correlated to GFR**  
Current methodology: Inulin, Iohexal, Creatinine Clearance, direct measurement
- **IDEXX SDMA was validated in rats** at the Reference Labs through a comparison to the gold standard method, liquid chromatography mass spectrometry.

## Sample submission


Assay	Species	Platform	Sample Type	Sample Volume*
SDMA	Rat	LCMS	Serum, Plasma	50 µL
SDMA	Mouse	LCMS	Serum, Plasma	50 µL
SDMA	NHP	LCMS	Serum, Plasma	100 µL
SDMA	Dog	LCMS	Serum, Plasma	100 µL
SDMA	Cat	Beckman	Serum, Plasma	100 µL
SDMA	Dog	Beckman	Serum, Plasma	100 µL
SDMA	Rat	Beckman	Serum, Plasma	100 µL

\* Minimum volume required


[See the SDMA Directory of Services for complete requirements](#)

# Kidneys in action

## Normal kidney function



## Kidney disease state



# Creatinine connection



## Does SDMA replace creatinine in diagnosing kidney disease?

SDMA does not replace creatinine in diagnosing kidney disease. While SDMA is a more reliable tool to evaluate kidney function/GFR and results should be considered first, creatinine is still complementary to SDMA for evaluating kidney function. A complete kidney evaluation should consist of a thorough history, physical examination, and evaluation of a minimum database, including CBC, chemistry profile with SDMA and electrolytes, and complete urinalysis.

## How does SDMA compare to creatinine when assessing kidney function?

As compared to creatinine, SDMA is a more reliable and sensitive indicator of kidney function in animals. SDMA increases earlier than creatinine in dogs and cats with kidney disease. SDMA increases with as little as 25% loss of kidney function (average 40% loss of kidney function). Creatinine does not increase until up to 75% of kidney function is lost. Creatinine is a breakdown product of muscle and is therefore impacted by lean body mass, whereas SDMA is not.

## GFR in relation to SDMA and Creatinine



Notice how many more SDMA are picking up change in SDMA at the RI cutoff for abnormal

GFR < 2 mL/min/kg this would be reduced or abnormal for canine species (3m-4ml/min/kg is normal GFR)



# LCMS — the gold standard

## What is the gold standard method for measuring SDMA?

The gold standard methodology for measuring SDMA is Liquid Chromatography Mass Spectrometry (LCMS), which has shown good correlation with glomerular filtration rate.

## What is liquid chromatography mass spectrometry?

Mass spectrometry, often referred to as mass spec, is an advanced analytical technique that is becoming increasingly popular in bioresearch. In a nutshell, mass spec is used to accurately identify and quantify molecules within a sample.

It is not often used for clinical work; however, it is increasingly relied upon for academic research and by regulatory authorities such as the Food and Drug Administration (FDA).

## For which species has SDMA been validated and can it be run on other species?

SDMA has been validated at the Reference Labs for dogs, cats, horses, and rats.

Projects validating SDMA and establishing reference intervals for other species are ongoing. At the Reference Lab, SDMA results will be provided on routine non-species-specific chemistry profiles.

## How was the SDMA reference interval established for rats at the Reference Labs?

IDEXX SDMA was validated in rats at the Reference Labs through a comparison to the gold standard method, liquid chromatography mass spectrometry. Following validation, a reference interval study on clinically healthy adult rats was performed.



# Human clinical studies

## Is SDMA used in human research studies?

SDMA has been studied in human beings and has been shown to correlate with GFR and to be a more sensitive marker for detecting early loss of kidney function than creatinine or eGFR (estimated GFR) calculations. At this time SDMA is not used routinely in human medicine because traditional methods of measuring SDMA with liquid chromatography-mass spectrometry (LCMS) are not very practical. SDMA currently appears to only be used during human clinical studies. **IDEXX is the only company globally to have successfully developed a high-throughput assay for SDMA and has been partnering with Yale University to help bring this technology to human medicine.**

- 18 studies N = 2,136 people
- SDMA correlated with inulin clearance
  - R = 0.85 (CI 0.76–0.91)
- SDMA correlated with creatinine
  - R = 0.75 (0.46–0.88)

- SDMA and Sepsis
- SDMA

Nephrol Dial Transplant (2006) 21: 2446–2451  
doi:10.1093/ndt/gfl292  
Advance Access publication 4 July 2006

## Original Article

### Symmetric dimethylarginine (SDMA) as endogenous marker of renal function—a meta-analysis

Jan T. Kielstein<sup>1,4</sup>, Shelley R. Salpeter<sup>2</sup>, Stefanie M. Bode-Boeger<sup>3</sup>, John P. Cooke<sup>1</sup> and Danilo Eliser<sup>4</sup>

<sup>1</sup>Department of Cardiovascular Medicine, <sup>2</sup>Department of Medicine, Stanford University School of Medicine, Stanford, CA, USA, <sup>3</sup>Institute for Clinical Pharmacology, Otto-von-Guericke University, Magdeburg and <sup>4</sup>Department of Internal Medicine, Division of Nephrology, Medical School Hannover, Germany

## Abstract

**Background.** Dosing of most drugs must be adapted in renal insufficiency, making accurate assessment of renal function essential in clinical medicine. Furthermore, even modest impairment of renal function has been recognized as a cardiovascular risk factor. The purpose of this analysis was to identify the role of symmetric dimethylarginine (SDMA), the structural isomer of the cardiovascular risk marker asymmetric dimethylarginine, as an endogenous marker of renal function.

**Methods.** Comprehensive searches of Medline and the Cochrane Library from 1970 to February 2006 were performed to identify studies that evaluated the correlation between SDMA and renal function. The search was augmented by scanning references of

## Symmetric dimethylarginine (SDMA)—a novel endogenous marker for renal dysfunction

The dosing of most drugs must be adapted in renal insufficiency, making accurate assessment of renal function a prerequisite in clinical medicine, especially in the elderly. Furthermore, even a modest decline in renal function has been recognized as a cardiovascular risk factor [1]. However, the diagnosis of kidney impairment is hampered by the lack of reliable markers of glomerular filtration rate (GFR). The gold standard, inulin clearance, is cumbersome and expensive, reducing its utility in clinical practice. This is also true for estimations using the clearance of radioisotopes. Therefore, in clinical practice,

Link: <https://academic.oup.com/ndt/article/21/9/2446/1939355>

- Foundational Study
- Background – two main forms of dimethylarginine – ADMA And SDMA
- ADMA – associated with inflammatory condition, not a good marker of renal function.
- SDMA : Continual production from cells, increased with decreasing renal function
- SDMA rises with reduction in GFR

## Asymmetrical dimethylarginine plasma clearance persists after acute total nephrectomy in rats

Katari A. Carello<sup>1</sup>, Steven E. Whitesell<sup>1</sup>, Mary C. Lloyd<sup>1</sup>, Scott S. Billecke<sup>1</sup> and Louis G. D'Alecy<sup>1,2</sup>

Departments of <sup>1</sup>Molecular and Integrative Physiology and <sup>2</sup>Surgery (Vascular), University of Michigan Medical School, Ann Arbor, Michigan; and <sup>3</sup>Department of Surgery, William Beaumont Hospital, Royal Oak, Michigan

[illegible]

ischemia may independently trigger the release of ADMA as suggested by indirect evidence for release of an endogenous NOS inhibitor in response to acute hindlimb ischemia (23–26). A recent review of Fliser et al. (19) traces the progression of ADMA, an innocent marker, to its current position as a leading candidate for a mediator of atherosclerotic complications in patients. In contrast, the physiological control, origins, and fates of this endogenous NOS inhibitor has yet to be fully defined.

In 1970, Kakimoto and Akazawa (27) first reported the mutual occurrence of ADMA and symmetrical DMA (SDMA) in urine and detected concentrations  $> 300 \mu\text{M}$  higher than that of arginine in urine. In contrast, the plasma ADMA and SDMA concentrations they reported were  $< 0.5\%$  (0.3  $\mu\text{M}$ ) of the concentration of arginine (55.8  $\mu\text{M}$ ) in the plasma, reflecting a dramatic urinary concentrating effect on ADMA and SDMA and unequivocally demonstrating their renal excretion.

By 1976, intravenous injections of radiolabeled, methylated arginine and detected concentrations  $> 300 \mu\text{M}$  in urine were catalogued in some way, whereas SDMA was largely unchanged in urine excreted (39). Ogawa et al. (47) later identified an enzyme that hydrolyzed ADMA to citrulline. This enzyme, dimethylarginine dimethylaminohydrolase (DDAH), has been identified in the heart, brain, liver, lung, and skeletal muscle and has been shown to have particularly high activity in the kidney (48). The enzyme has been purified from the kidney of SDMA in *Xiphophorus* set to emerge as the primary renal excretory

Link: <https://www.physiology.org/doi/full/10.1152/ajpheart.00208.2005>

## Mouse and Rat models

### SDMA present in renal tissue and in plasma/serum

Nephrectomized rats can be used as models for renal pathology

Pflügers Arch – Eur J Physiol (2000) 439:524–531  
Digital Object Identifier (DOI) 10.1007/s004249900220

## ORIGINAL ARTICLE

M. Al Banchaabouchi · B. Marescau · I. Possemiers  
R. D'Hooge · O. Levillain · P.P. De Deyn

### $N^G, N^G$ -Dimethylarginine and $N^G, N'^G$ -dimethylarginine in renal insufficiency

Received: 6 August 1999 / Received after revision: 19 October 1999 / Accepted: 8 November 1999 / Published online: 12 January 2000  
© Springer-Verlag 2000

**Abstract** Asymmetric  $N^G,N^G$ -dimethylarginine (ADMA) and symmetric  $N^G,N^G$ -dimethylarginine (SDMA) are basic endogenous amino acids with a guanidino group. Our renal distribution study of dimethylarginines clearly indicates that, in mouse and rat, ADMA and SDMA levels are most abundant as protein-incorporated com-

did not change significantly. Man, rat and mouse show similarities as well as differences in metabolism.

**Key words**  $N^G,N^G$ -Dimethylarginine (ADMA) ·  $N^G,N^G$ -dimethylarginine (SMDA) · Dimethylarginines · Renal failure · Species comparison

Link: <https://www.ncbi.nlm.nih.gov/pubmed/10764210>