

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT



Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)		
HCT (%)		
HGB (g/dL)		
MCV/VGM (fL)		
MCH/TCMH (pg)		
MCHC/CCMH (g/dL)		
PLT (K/μL)		
RDW-SD (fL)		
RDW-CV (%)		
PDW (fL)		
PCT (%)		
MPV (fL)		
P-LCR (%)		
WBC/GB (K/μL)		
NEUT# (K/μL)		
LYMPH# (K/μL)		
MONO# (K ⁹ /μL)		
EO# (K/μL)		
RBC-O (10 ⁶ /μL)		
PLT-O (K/μL)		
RETIC (K/μL)		
% RETIC (%)		

L2 (Normal / Normale / 正常濃度)

LOT

03080805



01/24/21

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)	4.40	4.13 - 4.67
HCT (%)	38.7	35.2 - 42.2
HGB (g/dL)	13.2	12.7 - 13.7
MCV/VGM (fL)	87.9	79.9 - 95.9
MCH/TCMH (pg)	30.0	26.3 - 33.7
MCHC/CCMH (g/dL)	34.1	28.0 - 40.2
PLT (K/μL)	234	189 - 279
RDW-SD (fL)	43.7	36.7 - 50.7
RDW-CV (%)	14.6	11.6 - 17.6
PDW (fL)	7.9	4.9 - 10.9
PCT (%)	0.20	0.12 - 0.28
MPV (fL)	8.6	5.6 - 11.6
P-LCR (%)	8.8	2.9 - 14.7
WBC/GB (K/μL)	6.93	6.18 - 7.68
NEUT# (K/μL)	3.29	2.23 - 4.35
LYMPH# (K/μL)	2.11	1.49 - 2.73
MONO# (K ⁹ /μL)	0.84	0.08 - 1.60
EO# (K/μL)	0.69	0.33 - 1.05
RBC-O (10 ⁶ /μL)	4.31	3.61 - 5.01
PLT-O (K/μL)	211	161 - 261
RETIC (K/μL)	122.8	92.8 - 152.8
% RETIC (%)	2.79	1.99 - 3.59