

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT 53280804  02/15/26

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁶ /μL)	2.18	2.03 - 2.33
HCT (%)	15.8	12.8 - 18.8
HGB (g/dL)	5.2	4.9 - 5.5
MCV/VGM (fL)	72.6	65.1 - 80.1
MCH/TCMH (pg)	23.9	20.1 - 27.7
MCHC/CCMH (g/dL)	32.9	21.7 - 44.1
PLT (K/μL)	54	24 - 84
RDW-SD (fL)	45.2	36.2 - 54.2
RDW-CV (%)	17.9	14.4 - 21.4
PDW (fL)	9.2	5.2 - 13.2
PCT (%)	0.05	0.03 - 0.07
MPV (fL)	8.7	5.7 - 11.7
P-LCR (%)	11.8	4.3 - 19.3
WBC/GB (K/μL)	3.05	2.65 - 3.45
NEUT# (K/μL)	1.37	0.96 - 1.78
LYMPH# (K/μL)	1.09	0.69 - 1.49
MONO# (K ⁹ /μL)	0.30	0.00 - 0.60
EO# (K/μL)	0.29	0.09 - 0.49
RBC-O (10 ⁶ /μL)	2.19	1.79 - 2.59
PLT-O (K/μL)	41	11 - 71
RETIC (K/μL)	148.4	108.4 - 188.4
% RETIC (%)	6.81	5.11 - 8.51

L2 (Normal / Normale / 正常濃度)

LOT 53280805  02/15/26

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁶ /μL)	4.48	4.21 - 4.75
HCT (%)	36.3	32.8 - 39.8
HGB (g/dL)	12.0	11.5 - 12.5
MCV/VGM (fL)	81.1	73.1 - 89.1
MCH/TCMH (pg)	26.8	23.1 - 30.5
MCHC/CCMH (g/dL)	33.1	27.0 - 39.2
PLT (K/μL)	211	166 - 256
RDW-SD (fL)	48.2	41.2 - 55.2
RDW-CV (%)	17.7	14.7 - 20.7
PDW (fL)	9.8	6.8 - 12.8
PCT (%)	0.20	0.12 - 0.28
MPV (fL)	9.6	6.6 - 12.6
P-LCR (%)	15.6	9.7 - 21.5
WBC/GB (K/μL)	7.22	6.47 - 7.97
NEUT# (K/μL)	3.53	2.47 - 4.59
LYMPH# (K/μL)	2.19	1.57 - 2.81
MONO# (K ⁹ /μL)	0.73	0.00 - 1.49
EO# (K/μL)	0.77	0.41 - 1.13
RBC-O (10 ⁶ /μL)	4.36	3.66 - 5.06
PLT-O (K/μL)	175	125 - 225
RETIC (K/μL)	126.9	96.9 - 156.9
% RETIC (%)	2.83	2.03 - 3.63