

- e-CHECK*** (XS) Quality Control for the **ProCyte Dx*** Analyzer
- e-CHECK*** (XS) contrôle qualité pour l'analyseur **ProCyte Dx***
- e-CHECK*** (XS) Control de calidad para el analizador **ProCyte Dx***
- e-CHECK*** (XS) Qualitätskontrolle für das **ProCyte Dx*** Analysegerät
- e-CHECK*** (XS) Controllo qualità per l'analizzatore **ProCyte Dx***
- e-CHECK*** (XS) IDEXXプロサイトDx*の精度管理

L1 (Low / Bas / Bajo / Niedrig / Basso / 低濃度)

LOT



Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)		
HCT (%)		
HGB (g/dL)		
MCV/VGM (fL)		
MCH/TCMH (pg)		
MCHC/CCMH (g/dL)		
PLT (K/μL)		
RDW-SD (fL)		
RDW-CV (%)		
PDW (fL)		
PCT (%)		
MPV (fL)		
P-LCR (%)		
WBC/GB (K/μL)		
NEUT# (K/μL)		
LYMPH# (K/μL)		
MONO# (K ⁹ /μL)		
EO# (K/μL)		
RBC-O (10 ⁶ /μL)		
PLT-O (K/μL)		
RETIC (K/μL)		
% RETIC (%)		

L2 (Normal / Normale / 正常濃度)

LOT

53290805



02/15/26

Parameter Paramètre Parámetro Parameter Parametro 項目	Target Valeur cible Objetivo Ergebnis Valore ideale 目標値	Expected Range Intervalle attendu Rango esperado Referenzbereich Intervallo atteso 基準範囲
RBC/GR (10 ⁹ /μL)	4.47	4.20 - 4.74
HCT (%)	36.2	32.7 - 39.7
HGB (g/dL)	12.1	11.6 - 12.6
MCV/VGM (fL)	80.9	72.9 - 88.9
MCH/TCMH (pg)	27.1	23.4 - 30.8
MCHC/CCMH (g/dL)	33.4	27.3 - 39.5
PLT (K/μL)	218	173 - 263
RDW-SD (fL)	48.0	41.0 - 55.0
RDW-CV (%)	17.8	14.8 - 20.8
PDW (fL)	9.8	6.8 - 12.8
PCT (%)	0.21	0.13 - 0.29
MPV (fL)	9.6	6.6 - 12.6
P-LCR (%)	15.3	9.4 - 21.2
WBC/GB (K/μL)	7.13	6.38 - 7.88
NEUT# (K/μL)	3.41	2.35 - 4.47
LYMPH# (K/μL)	2.23	1.61 - 2.85
MONO# (K ⁹ /μL)	0.75	0.00 - 1.51
EO# (K/μL)	0.74	0.38 - 1.10
RBC-O (10 ⁶ /μL)	4.38	3.68 - 5.08
PLT-O (K/μL)	176	126 - 226
RETIC (K/μL)	121.4	91.4 - 151.4
% RETIC (%)	2.72	1.92 - 3.52